

Certificate of Analysis

PRODUCT NAME 2 x Red PCR Master Mix
LOT NUMBER 100.805
CONCENTRATION 40 mU/μl Taq DNA Polymerase

FUNCTIONAL ASSAYS	SEPCIFICATION	RESULTS
4kb PCR	suitable	pass

PHYSICAL ASSAYS	SEPCIFICATION	RESULTS
Endonuclease Assay	no activity	no activity detected
Exonuclease Assay	no activity	no activity detected
Ribonuclease Assay	no activity	no activity detected

QC Acceptance Date: 01.06.2007



Expiry Date: 06/2008

Michael Thate, Quality Control